



Standard Screws, Locking Screws, Drill Bits & Taps

Screw Diameter										
Thread Diameter	1.5	2.0	2.4	2.7	3.5	4.5	4.0	2.4	2.7	3.5
Screw Type	Cortical						Cancellous	Locking Head Liberty Lock (Stainless)		
							Full Thread			
Drill Bit & Tap Diameter										
Drill Bit for Guiding Hole	1.5	2.0	2.4	2.7	3.5	4.5	4.0			
Drill Bit for Threaded Hole	1.1	1.5	1.8	2.0	2.5	3.2	2.5	1.8	2.0	2.8
Drive Type	1.5mm Hex		2.0mm Hex	2.5mm Hex		3.5mm Hex	2.5mm Hex	T8		T15